

ISTITUTO COMPRENSIVO

"G.LOMBARDO RADICE "

SIRACUSA



Funded by the
Erasmus+ Programme
of the European Union



FROM OLIVE OIL


A PRODUCT FROM OUR TRADITION

**HOMEMADE OIL
SOAP**



In the ancient times, the oil not suitable as food was not thrown away, but reused to prepare the soap, excellent both for the hygiene of the body and to clean and sanitize the laundry.

It is an excellent detergent, especially suitable for sensitive skin, as it does not contain aggressive agents, but only very few natural ingredients.



INGREDIENTS




500g olive oil



64g sodium
hydroxide



150g water




METHOD

There are two methods:



Warm Method



Cold Method



This method is the easiest and fastest

METHOD



Before you begin it's important to wear



GLOVES




**PROTECTIVE
GOGGLES**



MASK

INGREDIENTS



TOOLS

IMMERSION MIXER



SILICONE OR
PLASTIC
CONTAINER




METHOD


500 g Olive Oil



150 g Water




64 g Sodium Idroside




Weigh olive oil, water and caustic soda in plastic containers, reminding you always to wear gloves, goggles and a mask

METHOD

Pour into the container with the water,
the caustic soda (attention, never the
water in the soda, but the soda in the
water)



Stir slowly until the soda is completely
melted and the mixture will have become
transparent



METHOD

Add to the previous solution the olive oil




We use an immersion blender to quickly mix everything until the mixture gets thickened



METHOD

Pour the mixture into a mold lined with oven proof paper



Cover with a sheet of film and a cloth and leave to stand in a cool, dry place for 24 hours



METHOD

After removing from the container, cut the soap into smaller pieces and leave to stand for at least 40 days



OUR SOAP IS NOW READY FOR OUR USES

